

## HEARTSTART XLT DEFIBRILLATOR/ MONITOR



### THE DEFIBRILLATOR FOR PATIENTS IN TRANSPORT

Every caregiver knows that rapid time-to-defibrillation is critical when sudden cardiac arrest strikes. So Philips HeartStart XLT enables the first caregiver on-scene to deliver a lifesaving shock within seconds of arrival.

HeartStart XLT, a lightweight defibrillator/monitor for use by both EMTs and paramedics, switches from automated to manual operation with the press of a button, and charges to its highest energy setting within three seconds. In automated mode, clear, concise voice prompts guide users through the defibrillation process, while an easy-to-read LCD displays the ECG waveform. When switched to manual mode, HeartStart XLT is a manual defibrillator/monitor, offering true 1-2-3 operation. ACLS caregivers can select from a range of energy levels for defibrillation and

synchronized cardioversion. ECG monitoring provides key information needed for rapid intervention and thorough care. Optional non-invasive pacing and SpO<sub>2</sub> monitoring can be added.

HeartStart XLT features SMART Biphasic technology, Philips' patented waveform proven in multiple real-world studies to be superior to high-energy, monophasic waveforms while minimizing heart dysfunction.<sup>1, 2, 3</sup>

The XLT simplifies data management and reporting. All patient data, including continuous ECG and a variety of events, are stored on a data card. When downloaded to the Philips HeartStart Event Review data management system, data can be edited, stored, and printed as required for quality control and reporting.



### Easy to Use

- **Simple 1-2-3 Operation.** Makes defibrillation intuitive for all users.
- **AED Mode.** Voice and text prompts guide users through the defibrillation process.
- **3- or 5-Wire ECG Monitoring.** Activated in both AED and manual modes.
- **Multi-function Defibrillator Pads.** Feature clear illustrations for proper placement.
- **Backup AC Power.** AC power module for charging battery and operating the defibrillator.

### Smart Biphasic Advantage

- **Clinically Proven.** Our SMART Biphasic waveform technology uses low-energy that's highly effective in terminating VF.
- **Tailored Energy Delivery.** Measures chest impedance and delivers a low-energy shock based on individual needs.

### Lightweight, Compact, Durable

- **Compact.** Easily fits on a standard stretcher.
- **Rugged.** Withstands the rigors of pre-hospital use and patient transport.
- **Highly Visible Yellow Case.** Makes it easy to see in an emergency.

### A Proven Product from a Reliable Source

- **Training and Service.** Tailored to fit your needs.
- **Warranty Options.** One-year on-site service or five-year unit exchange (exchange available in North America only).
- **Dependability.** Part of the complete family of Philips resuscitation products used by healthcare providers for over 30 years.



## PHILIPS MEDICAL SYSTEMS IS PART OF ROYAL PHILIPS ELECTRONICS

On the web: [www.medical.philips.com](http://www.medical.philips.com)

Via e-mail: [medical@philips.com](mailto:medical@philips.com)

By fax: +31 40 27 64 887

#### By postal service:

Philips Medical Systems  
 Global Information Center  
 I.B.R.S. / C.C.R.I. Numéro 11088  
 5600 VC Eindhoven  
 Pays-Bas / The Netherlands

#### CORPORATE HEADQUARTERS

Philips Medical Systems  
 Nederland B.V.  
 Postbus 10.000  
 5680 DA Best  
 Netherlands

#### UNITED STATES

Andover, Massachusetts  
 (800) 934-7372

#### CANADA

Markham, Ontario  
 (800) 291-6743

#### EUROPE, MIDDLE EAST AND AFRICA

Böblingen, Germany  
 Fax: (49) 7031 463 1552

#### LATIN AMERICA HEADQUARTERS

Sunrise, Florida  
 (954) 835-2600

#### ASIA PACIFIC HEADQUARTERS

Wanchai, Hong Kong  
 (852) 2821 5888

© Philips Medical Systems

North America Corporation, 2002

All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Philips Medical Systems North America Corporation reserves the right to make changes in specifications or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in the Netherlands, may 2003  
 4522 982 92171

### REFERENCES

1. Bardy GH et al. "Multicenter Comparison of Truncated Biphasic Shocks and Standard Damped Sine Wave Monophasic Shocks for Transthoracic Ventricular Defibrillation". *Circulation* 1996;94:2507-2514.
2. Gliner BE, White RD. "Recurrence of Out-of-Hospital VF Following Low-Energy Biphasic and High-Energy Monophasic Defibrillation Shocks". *J American College of Cardiology* (Abstract) 1999;33:127A.
3. Reddy RK et al. "Transthoracic Ventricular Defibrillation Causes Fewer ECG ST-Segment Changes After Shock". *Ann Emerg Med* 1997;30:127-134.4.